



Bio Basic Inc.

A world leader in serving science

CERTIFICATE OF ANALYSIS

Product	AMP BUFFER CONCENTRATE
Grade	HIGH PURITY GRADE
Product Code	ADB0030
Formula	C ₄ H ₁₁ NO
MW	89.14
CAS#	124-68-5
Lot No	

Test Items	Specifications	Results
Appearance	White, slightly syrupy liquid	
Purity	>99.0%	
PH(0.1M,Water) @25C	11.0-12.0	

Storage: Store at room temperature. Protect from moisture.

QF 21 Rev 2015

20 Konrad Cres. Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com