



---

## CERTIFICATE OF ANALYSIS

---

|              |   |
|--------------|---|
| Product      | TE BUFFER (1x), pH8.0, STERILE SOLUTION   |
| Grade        | BIOTECHNOLOGY GRADE   |
| Product Code | D0123   |
| Mol. Formula | Not Applicable  |
| Mol. Weight  | Not Applicable  |
| Appearance   | Colorless, haze-free liquid   |
| Description  | A blend of Biotechnology Grade Tris Buffer and EDTA Disodium Salt, pH adjusted to 8.0 |
| Lot No       |   |

---

---

| Assay             | Specifications | Results |
|-------------------|----------------|---------|
| Conductivity@25°C | 775 -975 umhos |         |
| PH@25°C           | 7.90-8.10      |         |
| DNase, RNase      | None           |         |

---

Storage: Room Temperature.