



## CERTIFICATE OF ANALYSIS

---

Product	TotalBlot Nylon Membrane
Grade	N/A
Product Code	D0157-3
Formula	N/A
MW	N/A
CAS#	N/A
Lot No	

---

---

Test Items	Specifications	Results
Size	30cmX3m (1 Roll)	

---

Description: Positively charged nylon membranes for protein or nucleic acid blotting applications.

Storage: R.T.

Caution:

- Avoid contact with eyes, skin and clothing.
- Wash all areas of contact with copious amounts of water.
- See Material Safety Data Sheet for additional information.
- For laboratory or further manufacturing use only.
- Not intended for household use.

Laboratory Control