



CERTIFICATE OF ANALYSIS

Product	SSC BUFFER 20X LIQUID CONCENTRATE
Grade	BIOTECH GRADE
Product Code	D0623
Formula	N/A
MW	N/A
Lot No	

Test Items	Specifications	Results
Appearance	Clear and colorless liquid	
pH @25C	6.85-7.15	

Description: A pH adjusted buffer consisting of 3M Sodium Chloride & 0.3M Sodium Citrate

Storage: Store at room temperature. Protect from moisture.