



## CERTIFICATE OF ANALYSIS

---

Product	dNTP Mixture, (dATP, dCTP, dGTP, dUTP mixture)
Grade	Biotech
Formula	N/A
M.W.	N/A
CAS#	N/A
Product Code	DB2236
Description	A solution that is 2mM dATP, 2mM dCTP, 2mM dGTP, 4mM dUTP,
Batch No	

---

---

Test Items	Specifications	Results
Appearance	Clear and colorless solution	
Purity	>98.9% (HPLC)	
Dnase, Rnase	None detected	
pH	7.0 ± 0.1	

---

**Storage:**

Long term (infrequent use, 1-2 times per month): -70°C

Daily/Weekly use: -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes. These fluctuations can greatly alter product stability.