



CERTIFICATE OF ANALYSIS

Product	Peptone-C (peptonized milk)
Product Code	G292
Lot#	

Test Items	Specifications	Results
Amino Nitrogen (AN)	$\geq 1.9\%$	
Total Nitrogen(TN)	$\geq 6\%$	
Loss on Drying	$\leq 6\%$	
Ash(Residue on ignition)	$\leq 10\%$	
pH (2% Solution)	6.5-7.5	
<u>Growth Supporting Properties</u>		
Peptone Agar(USP)	Satisfactory	
Broth Formulation	Satisfactory	
<u>Standard Plate count</u>		
Molda and Yeast	<5000CFU/g	
Coliforms	<100CFU/g	
E.coli	Negative	
Salmonella	Negative	

Storage: R.T. Protect from moisture