



CERTIFICATE OF ANALYSIS

Product	L-Histidine, hydrochloride, monohydrate
Grade	FCC
Product Code	HE806
Formula	C ₆ H ₉ N ₃ O ₂ .HCl.H ₂ O
MW	209.63
CAS#	5934-29-2
Lot No	

Test Items	Specifications	Results
Appearance	White powder	
FCC	certified	
Identification	Pass	
Lead	≤0.0005%	
Loss on Drying	≤0.3%	
Purity (Dried)	98.5~101.5%	
Residue on Ignition	≤0.1%	
Specific Rotation (dried)	+8.5~+10.5°	

Storage: Store at room temperature. Protect from moisture.