



CERTIFICATE OF ANALYSIS

Product Name	PIPES, sodium salt
Grade	High purity
Product Code	PDB0434
Formula	$C_8H_{16.5}N_2 Na_{1.5}O_6S_2$
MW	335.34
CAS#	100037-69-2
Lot No	

Test Items	Specifications	Results
Appearance	White powder	
Assay (Titration, Dry basis)	≥99.0%	
Water Content	≤3.0%	
Infrared	Complies	
pH (1% Di H ₂ O, 25°C)	6.5 ~ 7.1	
Solubility (10% water)	Clear colorless solution	
A ₂₈₀ (10% water)	≤0.050	
A ₂₆₀ (10% water)	≤0.050	
Heavy metals (as Pb)	≤5ppm	

Storage: Store at room temperature. Protect from moisture