



CERTIFICATE OF ANALYSIS

Product	Pepstatin (Pepstatin A)
Grade	Ultra Pure
Product Code	PDJ694
Formula	C ₃₄ H ₆₃ N ₅ O ₉
MW	685.9
CAS#	26305-03-3
Lot No	

Test Items	Specifications	Results
Appearance	White powder	
Purity(HPLC)	≥98.0%	
Solubility	Soluble in acetic acid or methanol	

Suggested starting concentration: 0.7µg/ml (1µM)

Storage: 2~8°C. Protect from light.