



## CERTIFICATE OF ANALYSIS

---

Product	Sodium tungstate, dehydrate
Grade	High Purity
Product Code	S0898
Formula	$\text{Na}_2\text{WO}_4 \cdot 2\text{H}_2\text{O}$
MW	329.86
CAS#	10213-10-2
Lot No	

---

---

Test Items	Specifications	Results
Appearance	White fine crystalline powder	
Identification	Passes test	
Assay	99.0-101.0%	
Insoluble matter	$\leq 0.01\%$	
Titration free base	$\leq 0.02\%$	
Chloride (Cl)	$\leq 0.005\%$	
Nitrogen Compounds (as N)	$\leq 0.001\%$	
Molybdenum (Mo)	$\leq 0.001\%$	
Sulfate (SO <sub>4</sub> )	$\leq 0.01\%$	
Arsenic (As)	$\leq 5\text{ppm}$	
Heavy Metals and Iron (as Pb)	$\leq 0.001\%$	

---

Storage: 18-25°C