



## CERTIFICATE OF ANALYSIS

---

Product	YT BROTH AGAR POWDER
Product Code	SD7017
Appearance:	Light beige, free-flowing powder
Lot No	

---

---

Test Items	Specifications	Results
Conductivity (1:10)@25C	REPORT	
pH (6.6g/200ml water) @25C	6.6-7.0	
Solubility( 6.6g/200ml water)	PASS	
Supports Growth of Bacteria	PASS	

---

Component (g/L):

Yeast Extract	5.0g
Tryptone	8.0g
NaCl	5.0g
Agar	15.0g

Usage: 33.0g/L

Storage: Store at room temperature. Protect from moisture.