



## CERTIFICATE OF ANALYSIS

---

Product	L-Tartaric acid
Grade	Reagent
Product Code	TB0925
Formula	C <sub>4</sub> H <sub>6</sub> O <sub>6</sub>
MW	150.09
CAS#	87-69-4
Lot No	

---

---

Test Items	Specifications	Results
Appearance	White crystals or powders	
Assay	99.5 ~ 100.5%	
Specific rotation	+12.0° to +13.0°	
Heavy metals	≤10 ppm	
Residue on ignition	≤0.5%	
Loss on drying	≤0.5%	
Sulfate	≤150ppm	
Oxalate	≤350ppm	
Arsenic	≤2ppm	

---

Storage: R.T.