



CERTIFICATE OF ANALYSIS

Product	Thiourea
Grade	Reagent
Product Code	TB1943
Formula	(NH ₂) ₂ CS
MW	76.12
CAS#	62-56-6
Lot#	

Test Items	Specifications	Results
Appearance	White to Off-white solid	
Purity	≥99.0%	
Insoluble in water	≤0.002%	
Residue on ignition	≤0.005%	
Loss on drying	≤0.5%	
Thiocyanate	≤0.005%	
Transparency of solution	to pass test	

Storage: Store at room temperature.