



## CERTIFICATE OF ANALYSIS

---

|              |   |
|--------------|---|
| Product      | TT-SDS ( Tris-Tricine-SDS buffer) Premix powder |
| Grade        | Biotech   |
| Product Code | TD8135  |
| Formula      | N/A   |
| MW           | N/A   |
| CAS#         | N/A   |
| Lot No       |   |

---

---

| Test Items                        | Specifications                            | Results |
|-----------------------------------|---|---------|
| A single strength buffer contains | 0.01M Tris,<br>0.01M Tricine<br>0.01% SDS |         |

---

Each pack makes 1 liter of 10x TT-SDS buffer or 10 liters of 1XTT-SDS buffer.

Description: Tris-Tricine (TT)-SDS Buffer may replace Glycine with Tricine for use in standard TG and TG-SDS buffer system, has better resolution of low molecular weight proteins & peptides ( less than 10 kDa).

Storage: R.T.