



CERTIFICATE OF ANALYSIS

Product	Myokinase
Grade	N/A
Product Code	Y46622005
Formula	N/A
MW	N/A
CAS#	9013-02-9
Lot No	

Test Items	Specifications	Results
Source	from yeast	
Activity	Report	
Specific Activity	>200 units/mg protein	
Contaminants (as MK activity=100%)		
ATPase	< 0.01 %	
Phosphoglycerate kinase	< 0.1 %	
State		
Purified enzyme in solution		
One ml of solution contains the following		
Protein	1.0 ~ 30mg	
K-PO ₄	0.05 ~ 0.1 mol/L	
Glycerol	50%	
pH	7.0 ± 0.5	

Storage: 50% Glycerol solution at -20°C (do not freeze).