



Bio Basic Inc.

A world leader in serving science

Product Information

EZ-LINE Conical Centrifuge Tubes

Product information for SP61161:

Description

Conical centrifuge tubes with screw-caps. Tube exterior features volumetric gradation labels and a labeling area. Suitable for high-speed centrifuge work. Rack not included.

Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- External gradations indicate volume
- High strength propylene construction for high speed centrifugation



Item No.	Description	Package
SP61161	EZ-LINE Conical Centrifuge Tube, 50ml, 30x115mm, RCF 14,000 x g sterile, 25 per pack, 500 per case	500/case

Technical Specifications

Item No.	SP61161
Quantity per case	500/case
Packaging	25 Pack
RCF Rating	14,000 x g
Volume	50mL
Exterior Diameter	30mm
Length	115mm
Color	Clear Tube, Green Cap
Plastic materials	Polypropylene
Sterile	Yes

Storage

Room temperature.

OF 24 Rev 0

20 Konrad Cres. Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com