



CERTIFICATE OF ANALYSIS


Product	Tris HCl Buffer 1 M solution, Sterile pH 8.0
Grade	High Purity
Product Code	SD8127
Formula	N/A
MW	N/A
CAS#	N/A
Lot No	L807R0J
Manufacture Date	August, 2017
Retest Date	August, 2019

Test Items	Specifications	Results
pH @ 25°C	7.9-8.1	Conforms
Tris	0.98-1.02M	Conforms
DNase	Non-Detected	Conforms

This product does NOT contain any animal derived components.

Storage: 18~25°C.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature: 
Title: General Manager
Date: August 07, 2017

QF 21 Rev 2015