



Product Information

EZ-LINE Conical Centrifuge Tubes

Product information for SP61126:

Description

Conical centrifuge tubes with screw-caps. Tube exterior features volumetric gradation labels and a labeling area. Suitable for high-speed centrifuge work.

Features:

- Sterilized with gamma irradiation
- Non-pyrogenic
- External gradations indicate volume
- High strength propylene construction for high speed centrifugation



Item No.	Description	Package
SP61126	EZ-LINE Conical Centrifuge Tube, PP/HDPE, 15ml, 17x120mm, RCF 13000xg, Sterile to SAL 10-6, 50 per pack, 500 per case	500/case

Technical Specifications

Item No.	SP61126
Quantity per case	500/case
Packaging	50 Pack
Volume	15mL
Exterior Diameter	17mm
Length	120mm
Color	Clear Tube, Green Cap
Plastic materials	Polypropylene
Sterile	Yes

Storage

Room temperature.