



CERTIFICATE OF ANALYSIS

| | |
|--------------|--------------------------------|
| Product | EDTA Buffer, 0,5M , 1X, pH 8.0 |
| Grade | Cell culture |
| Product Code | USD8135 |
| Lot No | |

| Test Items | Specifications | Results |
|------------|----------------|---------|
| Appearance | Clear solution | |
| pH @ 25°C | 7.8-8.2 | |
| Sterility | Pass | |
| DNase | None Detected | |
| RNase | None Detected | |
| Protease | None Detected | |

Storage: 18~25°C