



Product Information Blotting Paper

Product information for PP3322/PP3323/PP3324:

Description

Blotting Paper made from pure raw materials and used as a membrane gel support by its high absorbency. This product is suitable for the lysis/denaturation of colony or plaque lifts and Western blots.

Features:

- 100% cotton fiber, white, glossy
- Manufactured and tested specifically for chromatographic and blotting techniques. This ensures the wicking capability and uniformity of capillary action that are important in obtaining clean and even transfers during blotting
- Convenient sizes available in sheets precisely cut to the most popular gel and transfer membrane sizes. Allows 'out-of-the-box' usage and eliminates sheet-to-sheet variations

Item No.	W X L	Thickness	Package
PP3322	10 x 15cm	180g/m ²	100 sheet
PP3323	15 x 15cm	180g/m ²	100 sheet
PP3324	20 x 20cm	180g/m ²	100 sheet

Storage

Room temperature.

Other Notes

Products are intended for basic scientific research only and **NOT** intended for diagnostic use or resale.