

## Product Information

### 15ml Conical Centrifuge Tube Graduated with Screw Cap

#### Product information for BT888 series:

##### Description:

15ml Polypropylene Conical Centrifuge Graduated Tubes with screw caps. Unlike polystyrene tube, Polypropylene tubes are chemical inert withstand centrifugation speed of 5000g. Autoclavable. Tubes with screw cap are suitable for cell culture.

##### Features:

- Polypropylene tube
- With screw cap
- Less than 5000g for centrifugation speed
- Suitable for cell culture
- Autoclavable
- Storage RT



##### Technical Specifications:

Cat.No.	BT888-N	BT888-SN	BT888-SB	BT888-SR	BT888-SY
Description	Tube & Natural Cap	Tube & Natural Cap	Tube & Blue Cap	Tube & Red Cap	Tube & Yellow Cap
Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Capacity(ml)	15	15	15	15	15
Bottom Style	Conical	Conical	Conical	Conical	Conical
Neck Style	Cylinder	Cylinder	Cylinder	Cylinder	Cylinder
Top Style	Threaded	Threaded	Threaded	Threaded	Threaded
Cap Style	Screw Cap	Screw Cap	Screw Cap	Screw Cap	Screw Cap
Self-Standing	No	No	No	No	No
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C
Graduations	Yes	Yes	Yes	Yes	Yes
Graduation Range	1.5~14.5ml	1.5~14.5ml	1.5~14.5ml	1.5~14.5ml	1.5~14.5ml
Graduation Interval	0.5ml	0.5ml	0.5ml	0.5ml	0.5ml
Sterile	No	Yes	Yes	Yes	Yes
Dnase-free	No	Yes	Yes	Yes	Yes
Rnase-free	No	Yes	Yes	Yes	Yes
Max RCF	5000g	5000g	5000g	5000g	5000g
Cap Diameter	21.67 mm	21.67 mm	21.67 mm	21.67 mm	21.67 mm
Outer Diameter	16.5 mm	16.5 mm	16.5 mm	16.5 mm	16.5 mm
Length	120 mm	120 mm	120 mm	120 mm	120 mm