



CERTIFICATE OF ANALYSIS

Product	Hygromycin B, 50mg/ml solution (in PBS), sterile
Grade	USP equivalent
Product Code	BS725
Formula	N/A
MW	N/A
CAS#	N/A
Lot No	

Test Items	Specifications	Results
Appearance	Liquid	
Hygromycin B	50mg/ml \pm 2%	
Activity(as dry basis or Raw material)	\geq 1200U/mg	
Purity of Raw material (HPLC)	\geq 92%	
Recommended retest period	12 months	

Storage: 4°C. Keep dry. Warm to room temperature before opening.