

## **Product Information**

### **50ml Conical with Screw Cap**

#### **Product information for CT788 series:**

#### **Description:**

50ml Polypropylene Conical with screw caps. Unlike polystyrene tube, Polypropylene tubes are chemical inert withstand centrifugation speed of 3000RPM. Autoclavable. Tubes with screw cap are suitable for cell culture.

#### **Features:**

- Polypropylene tube
- With screw cap
- Withstand centrifugation speed of 3000RPM
- Suitable for cell culture
- Autoclavable
- Storage RT



#### **Technical Specifications:**

Cat.No.	CT788-Y	CT788-G	CT788-YS	CT788-GS
Description	Tube & Yellow Cap	Tube & Green Cap	Tube & Yellow Cap	Tube & Green Cap
Material	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Capacity(ml)	50	50	50	50
Bottom Style	Conical	Conical	Conical	Conical
Neck Style	Cylinder	Cylinder	Cylinder	Cylinder
Top Style	Threaded	Threaded	Threaded	Threaded
Cap Style	Screw Cap	Screw Cap	Screw Cap	Screw Cap
Self-Standing	No	No	No	No
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C	-80°C~121°C
Graduations	Yes	Yes	Yes	Yes
Graduation Range	5~47.5ml	5~47.5ml	5~47.5ml	5~47.5ml
Graduation Interval	2.5ml	2.5ml	2.5ml	2.5ml
Sterile	No	No	Yes	Yes
Dnase-free	No	No	Yes	Yes
Rnase-free	No	No	Yes	Yes
Max RPM	3000RPM	3000RPM	3000RPM	3000RPM
Cap Diameter	33.3 mm	33.3 mm	33.3 mm	33.3 mm
Outer Diameter	29 mm	29 mm	29 mm	29 mm
Length	117 mm	117 mm	117 mm	117 mm