


## 12-well Dual System Electrophoresis Apparatus

(Cat.#EDK614)



Run the full spectrum of horizontal electrophoresis experiments with this versatile package! EDK614 supports one or two student groups in two standard length gel trays for experiments that require less separation, or one long gel tray for experiments that require more. Produces excellent results in 30-40 minutes. Includes (2) 7 X 10cm Trays, (1) 7 X 14cm Tray, (2) 6/8 tooth combs.

### Features:

- Large Color Coded Push Tabs for Easy Lid Insertion & Removal
- Pour Spout for Buffer Disposal
- Improved Ventilation Reduces Lid Condensation
- User Replaceable Electrodes



For 1-2 Lab Groups



Cat. #EDK614

