



## CERTIFICATE OF ANALYSIS

---

Product	Lactate dehydrogenase
Grade	Reagent
Product Code	L0253
Formula	N/A
MW	N/A
CAS#	N/A
Lot No	

---

---

Test Items	Specifications	Specifications
Appearance	Freeze Dried Powder	
Activity	≥100 U/mg solid	
Protein	Report	
Specific Activity	~350U/mg protein	
Solubility (1% W/V in water)	Soluble in distilled water or dilute buffer	

---

Storage: -15~ -20°C