



PRODUCT INFORMATION

Nutrient Agar (SD7041)

Product Name: Nutrient Agar

Product Code: SD7041

Ingredients:

Tryptone Equivalent = 5 g

Yeast Extract = 1.5 g

Beef Extract = 1.5 g

Sodium Chloride = 5 g

Agar = 15g

Storage:

Dehydrated media should be stored in a dry area at 18-26°C.

Prepared media should be stored at 2-8°C.

Method of Preparation:

1. Suspend 28g of powder into 1L of deionized water.
2. Allow the powder to dissolve. Heat and stir gently if required.
3. Autoclave at 121°C for 15 minutes.
4. Allow the solution to cool to 50°C or below before adding antibiotics or additional nutritional supplements.