

Product information

QF 24 V4
V3 Sept 2020

Specimen Collection Swab, Sterile

Catalog #: Z006553-0001 / Z006554-0001
Size: 1 UNIT
Storage: 15-30°C*

*: Check storage conditions.

Product Description:

Bio Basic offers CE and FDA certified swabs for the collection of nasal (Z006553-0001) and oral (Z006554-0001) samples.

Specifications:

Specimen Collection Swab, Sterile, Nasal (Z006553-0001)

Swab Length	Flocked Tip Length	Flocked Tip Diameter	Breakpoint
145mm±5mm	20mm	1.8mm±0.2mm	85mm±5mm



Z006553-0001: Specimen Collection Swab, Sterile, Nasal

Specimen Collection Swab, Sterile, Oral (Z006554-0001)

Swab Length	Flocked Tip Length	Flocked Tip Diameter	Breakpoint
150±2mm	22mm	2.8±0.2mm	78mm



Z006554-0001: Specimen Collection Swab, Sterile, Oral

Storage:

Store unused swabs dry between 15-30°C.