

## Product information

QF 24 V4  
V1 June 2020

# 96-well 2ml Deep-well Plate for Magnetic Separation Device

**Catalog #:** Z006670-0001  
**Size:** 1 UNIT  
**Storage:** 18 to 25°C

\*: Product shipped and stored at room temperature.

### Features:

- PP material.
- Dimension: 127.76mm x 85.48mm x 45.60mm.

Item No.	Color	Sterile or Not	Package
Z006670-0001	Natural	Sterile	Individual

### Storage:

Room temperature.

