

Product information

QF 24 V4
 V1 July 2021

Digital Water Baths

Catalog #: BCM1400/1401/1402/1403
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Digital water baths, with gabled lid. Gabled lid diverts condensation to front and back of bath. Lid is removable. BCM1400 has a flat, non-removable lid. 8L and 12L models have a drain for easy removal of water. All models have seamless, stainless steel interiors which are rust-resistant and have no corners to collect debris or contaminants. 2 year warranty.



Features:

- Hinged, removable lid (except for BCM1400 – lid cannot be removed)
- Seamless stainless steel chamber
- Digital controls

Ordering Information:

BCM1400 - 2L digital water bath, 115V
BCM1401 - 4L digital water bath, 115V
BCM1402 - 8L digital water bath, 115V
BCM1403 - 12L digital water bath, 115V

BCM1402 - 8L digital water bath, 115V

Technical Data:

	BCM1400	BCM1401	BCM1402	BCM1403
Temperature Range	Ambient +5°C to 100°C			
Temperature Accuracy	±0.2°C at 37 °C		±0.5°C at 37 °C	
Temperature Uniformity	±0.2°C at 37 °C			
Exterior Dimensions	6.5(W) x 8(D) x 6(H)''	7.7(W) x 13.9(D) x 9(H)''	17.5(W) x 10.5(D) x 12.4(H)''	18(W) x 10.5(D) x 16.2(H)''
Chamber Dimensions	5.5(W) x 6(D) x 4(H)''	5.3(W) x 9.4(D) x 5.5(H)''	12.8(W) x 10.4(D) x 5.9(H)''	12.8(W) x 12.4(D) x 7.9(H)''
Weight	4.5 lbs / 2kg	15 lbs / 6.5 kg	17 lbs / 8 kg	17 lbs / 8 kg
Electrical	115V, 50-60 Hz			
Warranty	2 years			