

Product information

QF 24 V4
 V1 July 2021

Infrared MicroSterilizer

Catalog #: BCM1543
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Instantly sterilizes various metal and borosilicate glass instruments. Ceramic heating element is protected by a stainless steel guard chamber. Minimum heat setting can be used as a standby mode to extend the MicroSterilizer's usable life. 1 year warranty.

Features:

- Sterilizes instruments in 5-7 seconds
- Ideal for inoculating loops and needles and glass tube/pipette mouths



Technical Data:

Max. Temperature	815°C/1500°F
Min. Temperature	500°C/900°F
Sterilization Time	5-7 seconds
Dimensions (WxDxH)	4 x 6 x 6.5 inches 10.2 x 15.2 x 16.5 centimeters
Weight	4lbs/1.3kg
Chamber Diameter	0.55 inches / 14mm
Chamber Length	6 inches / 153 mm
Electrical	115V, 50-60Hz, 170W
Warranty	1 Year



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.