

Product information

QF 24 V4
 V1 July 2021

Mini Vortexer

Catalog #: BCM1600
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Compact, powerful vortexer for even the most crowded lab bench. Unique spill-proof head keeps motor dry. Activates when pressure is applied to the blue head. 2 year warranty.

Features:

- Powerful motor
- Compact, palm-sized footprint
- Pressure-sensitive head instantly starts and stops when a sample is touched to the surface



Technical Data:

Speed	2,800rpm
Orbit	4mm
Weight	1lbs/0.4kg
Dimensions (WxDxH)	3.7 x 3.9 x 2.6 inches 9.4 x 9.9 x 6.6 centimeters
Electrical	Universal Voltage 100-240V, 50-60Hz, 0.7A
Activation	Touch
Operating Temperature	+4°C to +65°C
Warranty	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.