

## Product information

QF 24 V4  
 V1 July 2021

## Vortexers

**Catalog #:** BCM1601 / BCM1607  
**Size:** 1 UNIT  
**Storage:** 18-25°C

### Product Description:

Variable speed vortexers. BCM1601 includes a cup head. BCM1607 includes a multipurpose head. Optional heads are also available. Please see products BCM1602 through BCM1606. 2 year warranty.

### Features:

- Quiet motor
- Multipurpose head allows for standard vortexing, and holds up to 8 microtubes
- Spill-proof controls
- Cold room and incubator-compatible
- Additional heads available



**BCM1601 - Mixer Vortexer 115V -  
 Standard cup head**

### Technical Data:

<b>Speed Range</b>	200-3200rpm
<b>Operating Mode</b>	Touch or continuous
<b>Electrical</b>	115V, 50Hz, 180W
<b>Dimensions (WxDxH)</b>	5 x 6.3 x 6.7 inches 13 x 16 x 17 centimeters
<b>Orbit</b>	3mm
<b>Weight</b>	8.4lbs/3.8kg
<b>Operating Temp Range</b>	+4°C - +65°C
<b>Warranty</b>	2 Years

### Ordering Information:

**BCM1601** - Mixer Vortexer 115V - Standard cup head

**BCM1607** - Vortexer 115V – Multipurpose head



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.  
 NOT INTENDED FOR HUMAN OR ANIMAL USE.