

## Product information

QF 24 V4  
 V1 July 2021

# Tube Rocker

**Catalog #:** BCM1610  
**Size:** 1 UNIT  
**Storage:** 18-25°C

### Product Description:

Fixed speed tube rocker with grooved mat. 2 year warranty.

### Features:

- No tools or mounting hardware required
- Ideal for tubes 1.5ml to 50ml



**BCM1610 - Compact Rocker with Grooved Mat for Tubes, 115V**

### Technical Data:

<b>Speed</b>	24rpm
<b>Tilt Angle</b>	±30°
<b>Load Capacity</b>	5 lbs / 2.3 kg
<b>Dimensions (WxDxH)</b>	13.3 x 5.5 x 4 inches 33.6 x 13.5 x 10 centimeters
<b>Weight</b>	4 lbs / 1.8 kg
<b>Electrical</b>	115V, 50-60Hz, 50W
<b>Operating Temp Range</b>	+4°C to +65°C
<b>Platform Size</b>	12 x 4.5 inches / 30 x 11.4 cm
<b>Warranty</b>	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.  
 NOT INTENDED FOR HUMAN OR ANIMAL USE.