

Product information

QF 24 V4
 V1 July 2021

Platform Rocker

Catalog #: BCM1613
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Fixed speed platform rocker with flat mat. 2 year warranty.

Features:

- Ideal for blot washing and gel staining



BCM1613 - Platform Rocker, fixed speed,
 12" x 12" with flat mat, 115V

Technical Data:

Speed	12rpm
Tilt Angle	±8°
Load Capacity	5 lbs / 2.3 kg
Dimensions (WxDxH)	13.5 x 12 x 4.7 inches 34 x 30 x 12 centimeters
Weight	6.5lbs / 3kg
Electrical	115V, 50-60Hz, 50W
Operating Temp Range	+4°C to +65°C
Platform Size	12 x 12 inches / 30 x 30 cm
Warranty	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.