

Product information

QF 24 V4
 V1 July 2021

3D Nutating Shakers

Catalog #: BCM1615/1616
Size: 1 UNIT
Storage: 18-25°C

Product Description:

Multifunctional mixers that provide a 3D nutating motion for thorough mixing. Brushless, maintenance-free motors. 2 year warranty.

Features:

- Includes dimpled mat for use with centrifuge or blood tubes
- Stable tri-directional motion for uniform mixing
- Mat is autoclavable
- Compatible with BCM1645 bungee cords to secure trays to platform
- Safe for incubator or cold room use

Ordering Information:



BCM1616 - Mixer 3D Nutating Shaker, fixed speed, 12" x 12" platform, dimpled mat, 115V

BCM1615 - Mixer 3D Nutating Shaker, fixed speed, 10" x 7.5" platform, dimpled mat, 115V

BCM1616 - Mixer 3D Nutating Shaker, fixed speed, 12" x 12" platform, dimpled mat, 115V

Technical Data:

Speed	24rpm
Tilt Angle	±20°
Load Capacity	4 lbs / 1.8 kg
Weight	4.4 lbs/2 kg
Electrical	115V, 50-60Hz, 50W
Dimensions (WxDxH)	BCM1615: 10.5 x 7.5 x 8 inches 27 x 19 x 20 centimeters BCM1616: 12 x 12 x 8 inches 30 x 30 x 20 centimeters
Platform Size	BCM1615: 10.5 x 7.5 inches 27 x 19 centimeters BCM1616: 12 x 12 inches 30 x 30 centimeters
Operating Temperature	+4°C to +65°C
Warranty	2 Years



PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC RESEARCH ONLY.
 NOT INTENDED FOR HUMAN OR ANIMAL USE.