



## CERTIFICATE OF ANALYSIS

---

Product	Phenol, saturated (pH4.5)
Grade	High Purity Grade
Product Code	PC6920
Formula	N/A
MW	N/A
CAS#	108-95-2
Lot No	

---

---

Test Items	Specifications	Results
pH	4.5±0.2	
Saturated with	Tris-HCl	
Absorbance 330nm	≤ 0.2	
Absorbance 405nm	≤ 0.1	
Absorbance 510nm	≤ 0.1	

---

Storage: -15~-20°C