



CERTIFICATE OF ANALYSIS

Product	RNA-Be-Locked (Tissue) Reagent
Product Code	RT4171
LOT#	O102DR0J
Expiration date	July 01, 2021

Kit components


Components	RT4171, 25ml
RNA-Be-Locker A	25 ml
Protocol	1

This is to certify that units of the above mentioned lot number were tested and found to comply with the specifications of the grade listed. Bio Basic expressly disclaims all warranties, expressed or implied, including the implied warranties of merchantability and fitness for a particular purpose. The following are the actual analytical results obtained:

Result name	Units	Specification	Test Value
Appearance		Report	Kit
Function Test	Pass/Fail	=Pass Test	Pass test

Storage: Transportation and storage at room temperature. Valid for 2 year.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature: 
Title: General Manager
Date: January 02, 2020