



## CERTIFICATE OF ANALYSIS

Product	TBS buffer (Tris-Buffered Saline), 20x Solution
Product Code	UA0045
Grade	Cell culture

Test Items	Specifications	Results
Appearance	Clear Solution	Pass
pH (1X water at 25°C)	7.4 ±0.2	Pass
Sterile by 0.2 µM filtration	Comply	Pass
DNase	None Detected	Pass
RNase	None Detected	Pass
Protease	None Detected	Pass
Tris	0.5M	Pass
NaCl	3.0M	Pass

Storage: 18~25°C.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature:   
Title: General Manager  
Date: June 03, 2022