



## CERTIFICATE OF ANALYSIS

Product Code: AS88-60158-500mL

Product Name: Sodium phosphate, 0.5M buffer soln., pH 7.6

Lot No. :

Appearance: Clear colorless liquid

pH:  $7.6 \pm 0.15$

Sodium Phosphate dibasic/  
Sodium phosphate monobasic 0.5M

Method: All components are dissolved in 18.2 megohm-cm  $\pm 1$  water (CDI, 18 megohm, UV treated), followed by pH adjustment, and filtered through 0.22 micron filter.