



## CERTIFICATE OF ANALYSIS

Product Code: AS88-63682-1L

Product Name: TRIS-buffered saline (TBS, 20X), with 0.5% Tween 20

Lot No. : X23I500

Appearance:	Clear, viscous liquid
pH:	7.4 ± 0.15
Tris-HCl	1 M
NaCl	3 M
Tween-20	1 %

Method: Prepared in 18.2 megohms-cm ± 1 water and autoclaved 15 psi (121°C) for 30 minutes on liquid cycle.