



CERTIFICATE OF ANALYSIS

Product Code: AS88-67760-500mL

Product Name: Ringer's Solution, sodium free, lactate-buffered

Lot No. :

Appearance:	Colorless clear liquid
pH:	7.4 ± 0.15
Choline Chloride:	100 mM
Potassium Chloride:	3.7 mM
Magnesium Sulfate:	1.2 mM
Calcium Chloride:	1.3 mM
Choline bicarbonate:	28 mM
Lactic Acid:	20 mM

Method: Prepared in 18.2 megohms-cm ± 1 water and then filtered through a 0.22-micron filter.