



## CERTIFICATE OF ANALYSIS

Product Code: AS88-60473-500mL

Product Name: BSA blocking buffer, 3% in PBS, with azide

Appearance: Clear liquid

pH:  $7.4 \pm 0.15$

Sodium chloride: 137 mM

BSA in PBS: 3 %

KCl: 2.7 mM

Sodium azide: 0.02 %

Potassium Phosphate Monobasic: 2 mM

Sodium Phosphate Dibasic: 10 mM

Method: Prepared in  $18.2 \text{ megohm-cm} \pm 1$  low endotoxin water, pH adjusted to 7.4 and then filtered through 0.22-micron filter.