



## CERTIFICATE OF ANALYSIS

Product Code: AS88-60565-1L

Product Name: Sodium phosphate, 0.2M buffer soln., pH 9.0

Lot No. :

Appearance: Clear liquid

pH:  $9.0 \pm 0.15$

NaH<sub>2</sub>PO<sub>4</sub>/Na<sub>2</sub>HPO<sub>4</sub>: 200 mM

Method: Ultrapure sodium phosphates in given concentration was dissolved in 18.2 megohm-cm  $\pm$  1 water (CDI, 18 megohm, UV treated), adjusted to pH 9.0 and then filtered through 0.22- micron filter.