



CERTIFICATE OF ANALYSIS

Product Code: AS88-61089-250mL

Product Name: BSA blocking buffer, 5% in PBS

Lot No. :

Appearance: Clear viscous liquid

pH: 7.4 ± 0.15

Bovine Serum albumin: 5 %

Sodium Phosphate Dibasic: 10.1 mM

Potassium Phosphate Monobasic: 1.76 mM

Potassium Chloride: 2.7 mM

Sodium Chloride: 137 mM

Method: Prepared in $18.2 \text{ megohm-cm} \pm 1$ water and filtered through 0.22-micron filter.