



## CERTIFICATE OF ANALYSIS

Product Code: AS88-61288-250mL

Product Name: Sodium acetate, 3M aq. soln., pH 4.5, autoclaved

Lot No. :

Appearance: Clear liquid

pH:  $4.5 \pm 0.15$

Sodium acetate: 3 M

Method: All components were dissolved in 18.2 megohm-cm  $\pm 1$  water and autoclaved at 15 psi (121°C) for 30 minutes on liquid cycle and/or filtered through 0.22-micron filter.