



## CERTIFICATE OF ANALYSIS

Product Code: AS88-61928-250mL

Product Name: Sodium acetate, 3M aq. soln., pH 5.2, RNase free

Lot No. :

Appearance: Clear liquid

pH:  $5.2 \pm 0.15$

Sodium Acetate: 3 M

DEPC: 0.01 %

Method: Prepared in 18.2 megohm-cm  $\pm 1$  DNase/RNase free water and autoclaved at 15 psi (121°C) for 30 minutes on liquid cycle and/or filtered through 0.22-micron filter.