



CERTIFICATE OF ANALYSIS

Product Code: AS88-61934-100mL

Product Name: Sodium acetate, 3M aq. soln., pH 7.0, RNase free

Lot No. :

Appearance: Clear liquid

pH: 7.0 ± 0.15

Sodium Acetate: 3 M

DEPC: 0.01 %

Method: Prepared in 18.2 megohm-cm ± 1 DNase/RNase free water and autoclaved at 15 psi (121°C) for 30 minutes on liquid cycle and/or filtered through 0.22-micron filter.