



CERTIFICATE OF ANALYSIS

Product Code: AS88-62305-500mL

Product Name: BSA blocking buffer, 5% in TBS, with 0.02% sodium azide

Lot No. :

Appearance: Clear viscous liquid

pH: 7.4 ± 0.15

BSA: 5 %

Sodium Azide: 0.02 %

TRIS: 50 mM

Sodium Chloride: 137mM

Method: Prepared in $18.2 \text{ megohms-cm} \pm 1$ and then filtered through a 0.22-micron filter.