



## CERTIFICATE OF ANALYSIS

Product Code: AS88-63283-QA

Product Name: BCA protein assay kit

Lot No. :

Appearance: Contains clear and colored liquids

Solution A

Sodium bicinchoninate: 1 % (w/v)

Sodium carbonate: 2 % (w/v)

Sodium bicarbonate: 0.95 % (w/v)

Sodium tartrate: 0.16 % (w/v)

(dibasic dehydrate)

Sodium hydroxide: 0.4 % (w/v)

Solution B

Cupric sulfate (5X hydrated): 4 % (w/v)

Method: Prepared in 18.2 megohm-cm  $\pm$  1 water and filtered through 1-micron filter.