




CERTIFICATE OF ANALYSIS

Product	Pepsin, 1:3000
Grade	USP Grade
Product Code	PB0411
Source	Porcine Stomach
CAS#	9001-75-6
Lot No	

Test Items	Specifications	Results
Description	White or light yellow powder, no impurities, with smell of enzyme application	
Enzyme activity	$\geq 3,700$ NF.units/mg	
Lead (Pb)	≤ 5 mg/kg	
Total arsenic	≤ 3 mg/kg	
Total number of colonies	≤ 50000 cfu/g	
Coliform	≤ 30 cfu/g	
Escherichia coli	≤ 3 mpn/g	
Salmonella (25g)	Not detected	
Net weight	Complies	

Storage: 2-8°C. Preserve in a well-closed container and keep in dry, dark place.

This Signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature: 
Title: General Manager
Date: October 25, 2023